

Masimo® Product Catalog


2012



Masimo Corporation
40 Parker
Irvine, CA 92618
Tel: (949) 297-7000
Fax: (949) 297-7001
www.masimo.com

© 2012 by Masimo Corporation.

This catalog contains Masimo proprietary material and is provided for information purposes only. No part of this catalog may be reproduced or retransmitted without Masimo's prior written consent. Product availability is subject to change without notice.







Masimo, Rainbow, SpCO[®], SpHb[™], SpMet[®], Cleanshield, DCI, LNCS, LNOP, Rad, Radical, , SatShare, SET, Signal Extraction Technology, Signal I.Q., Rad-5, Rad-8, Rad-9, RadLink are registered trademarks of Masimo Corporation.





M-LNCS, Blue, Radical-7, Rad-57, Pulse CO-Oximeter, Red, TF-I, RadNet, TrendCom and PPO+ are trademarks of Masimo Corporation.




The devices advertised in the catalog may not have been licensed in accordance with applicable Federal Laws.



CONFIDENTIAL AND PROPRIETARY

Monitoring Devices


Masimo SET® Pulse Oximetry Devices				
Masimo PN	Description	UM	Product Image	Notes
9021	Radical® Docking Station, RDS-1 <i>Comes with a choice of power cord and serial output connector, Analog output / nurse call connector and SatShare cable connector.</i>	EA		
1311	Radical® Docking Station, RDS-1B <i>With Extended Battery comes with a choice of power cord and serial output connector, Analog output / nurse call connector and SatShare cable connector.</i>	EA		
1392	Radical® Docking Station, RDS-2 <i>Comes with choice of Power cord (International Only), and no communication capabilities.</i>	EA		
9023	Radical® Docking Station, RDS-3 limited communication capabilities <i>Comes with choice of power cord and serial output connector and Analog output / nurse call connector.</i>	EA		
9190	Rad-8®, Horizontal <i>Comes with either an LNOP or LNCS Patient Cable and a power cord.</i>	EA		
9191	Rad-8®, Vertical <i>Comes with either an LNOP or LNCS Patient Cable and a power cord.</i>	EA		

Masimo SET® Pulse Oximetry Devices				
Masimo PN	Description	UM	Product Image	Notes
9230	Radical® (20-pin) Handheld, Color Screen <i>Includes SpO₂, Pulse Rate and Perfusion Index (PI) parameters as well as one Red LNOP, Red LNCS or Rainbow Patient Cable.</i>	EA		Does not include Rainbow Activation License. Device cannot be upgraded with Rainbow parameters.
9196	Rad-5® Pulse Oximeter <i>Comes with a Red DCI-DC3 Sensor and choice of patient cable, a choice of available Colors of Handheld Boots and a choice of user manual.</i>	EA	 Shown in optional protective boot and with DCSC Sensor	Please refer to Handheld Device accessories for available colors.
9198	Rad-5® Pulse Oximeter with Sample Pack <i>Comes with choice of LNCS or LNOP Sensor Sample Pack and LNOP Patient Cable, a choice of available Colors of Handheld Boots and a choice of user manual.</i>			
9214	Rad-5® Pulse Oximeter (International Only) <i>Comes with an LNOP DCSC Sensor, a choice of available Colors of Handheld Boots and a choice of user manual.</i>			
9197	Rad-5v™ Pulse Oximeter <i>Comes with a Red DCI-DC3 Sensor and choice of patient cable, a choice of available Colors of Handheld Boots and a choice of user manual.</i>	EA	 Shown in optional protective boot and with DCSC Sensor	Please refer to Handheld Device accessories for available colors.
9199	Rad-5v™ Pulse Oximeter with Sample Pack <i>Comes with choice of LNCS or LNOP Sensor Sample Pack and LNOP Patient Cable, a choice of available Colors of Handheld Boots and a choice of user manual.</i>			
9215	Rad-5v™ Pulse Oximeter (International Only) <i>Comes with an LNOP DCSC Sensor, a choice of available Colors of Handheld Boots and a choice of user manual.</i>			
9500	Radical-7™ Touch Screen Handheld <i>Includes SpO₂, Pulse Rate and Perfusion Index (PI) parameters as well as one Red LNOP, Red LNCS or Rainbow Patient Cable. Can be configured with optional Rainbow parameters and features: Total Hemoglobin (SpHb™), Oxygen Content (SpOC™), Methemoglobin (SpMet®), Carboxyhemoglobin (SpCO®), Pleth Variability Index (PVI™). This unit ships standard with a 802.11 wireless radio.</i> <i>Can be loaded with Respiration Rate from Pleth (RRp) and Adaptive Threshold Alarms (ATA™) in cleared countries only. Please check with customer service on availability.</i>	EA		Optional Rainbow parameters and features are sold separately and can be added to this unit upon request. Please refer to Rainbow SET Device options for parameter/feature part numbers. The Radical-7 is compatible with all Radical Docking Stations.


Masimo Rainbow SET® Devices				
Masimo PN	Description	UM	Product Image	Notes
9086	Radical-7™ Handheld, Color Screen Includes SpO ₂ , Pulse Rate and Perfusion Index (PI) parameters as well as one Red LNOP, Red LNCS or Rainbow Patient Cable. Can be configured with optional Rainbow parameters and features: Total Hemoglobin (SpHb™), Oxygen Content (SpOC™), Methemoglobin (SpMet®), Carboxyhemoglobin (SpCO®), Pleth Variability Index (PVI™), Desat Index 3D alarm and PI Delta 3D alarm.	EA		Optional Rainbow parameters and features are sold separately and can be added to this unit upon request. Please refer to Rainbow SET Device options for parameter/feature part numbers. The Radical-7 is compatible with all Radical Docking Stations.
9153	Radical-7™ Handheld, Blue Screen Includes SpO ₂ , Pulse Rate and Perfusion Index (PI) parameters as well as one Red LNOP, Red LNCS or Rainbow Patient Cable. Can be configured with optional Rainbow parameters and features: Total Hemoglobin (SpHb), Oxygen Content (SpOC), Methemoglobin (SpMet), Carboxyhemoglobin (SpCO), Pleth Variability Index (PVI), Desat Index 3D alarm and PI Delta 3D alarm.	EA		
9226	Rad-87™, Horizontal Includes SpO ₂ , Pulse Rate and Perfusion Index (PI) parameters. Includes either Red LNCS, Red LNOP or Rainbow RC Patient Cable. Can be configured with optional Rainbow parameters and features: Total Hemoglobin (SpHb), Oxygen Content (SpOC), Methemoglobin (SpMet), Carboxyhemoglobin (SpCO), Pleth Variability Index (PVI), Acoustic Respiration (RRa), Adaptive Threshold Alarm (ATA), in countries where above parameters are cleared. Comes with serial output connector and Analog output / Nurse call connector.	EA		Rainbow Acoustic Monitoring capable. Optional Rainbow parameters and features are sold separately and can be added to this unit upon request. Please refer to Rainbow SET Device options for parameter/feature part numbers.
9227	Rad-87™, Horizontal with Wireless Radio Includes SpO ₂ , Pulse Rate and Perfusion Index (PI) parameters and integrated wireless radio. Includes either Red LNCS, Red LNOP or Rainbow RC Patient Cable. Can be configured with optional Rainbow parameters and features: Total Hemoglobin (SpHb), Oxygen Content (SpOC), Methemoglobin (SpMet), Carboxyhemoglobin (SpCO), Pleth Variability Index (PVI), Acoustic Respiration (RRa), Adaptive Threshold Alarm (ATA), in countries where above parameters are cleared. Comes with serial output connector and Analog output / Nurse call connector.			

Masimo Rainbow SET® Devices				
Masimo PN	Description	UM	Product Image	Notes
9228	Rad-87™, Vertical Includes SpO ₂ , Pulse Rate and Perfusion Index (PI) parameters. Includes either a Red LNCS, Red LNOP or Rainbow RC Patient Cable. Can be configured with optional Rainbow parameters and features: Total Hemoglobin (SpHb), Oxygen Content (SpOC), Methemoglobin (SpMet), Carboxyhemoglobin (SpCO), Pleth Variability Index (PVI), Acoustic Respiration (RRa), Adaptive Threshold Alarm (ATA), in countries where above parameters are cleared. Comes with serial output connector and Analog output / Nurse call connector.	EA		Rainbow Acoustic Monitoring capable. Optional Rainbow parameters and features are sold separately and can be added to this unit upon request. Please refer to Rainbow SET Device options for parameter/feature part numbers.
9229	Rad-87™, Vertical with Wireless Radio Includes SpO ₂ , Pulse Rate and Perfusion Index (PI) parameters and integrated wireless radio. Includes either a Red LNCS, Red LNOP or Rainbow RC Patient Cable. Can be configured with optional Rainbow parameters and features: Total Hemoglobin (SpHb), Oxygen Content (SpOC), Methemoglobin (SpMet), Carboxyhemoglobin (SpCO), Pleth Variability Index (PVI), Acoustic Respiration (RRa), Adaptive Threshold Alarm (ATA), in countries where above parameters are cleared. Comes with serial output connector and Analog output / Nurse call connector.			
9221	Pronto-7™ Spot Check Pulse CO-Oximeter Includes SpO ₂ , Pulse Rate, Perfusion Index (PI) and SpHb parameters. Rainbow 4D DC Sensor not included.	EA		Rainbow 4D DC Reusable Sensor must be ordered separately.
9167	Pronto Spot Check Pulse CO-Oximeter Includes SpO ₂ , Pulse Rate, Perfusion Index (PI) and SpHb parameters. Rainbow Sensor not included.	EA		Rainbow Sensor must be ordered separately.

Masimo Rainbow SET® Devices

Masimo PN	Description	UM	Product Image	Notes
9216 - U 9216	<p>Rad-57™ Configurable Pulse CO-Oximeter</p> <p>Includes SpO₂, Pulse Rate and Perfusion Index (PI) parameters.</p> <p>The Rad-57 can ship with or be field upgraded with Carboxyhemoglobin (SpCO), Methemoglobin (SpMet), Total Hemoglobin (SpHb), Oxygen Content (SpOC), and/or PVI</p> <p>When shipped with Rainbow parameters, Rainbow Sensor must be purchased separately.</p> <p>Rainbow Sensors not included.</p>	EA	 <p>Shown with Rainbow Adult Reusable Sensor and in optional protective boot</p>	<p>Please refer to Handheld Device accessories for available colors.</p> <p>9216-U is for SpO₂ only Devices.</p>




Masimo Rainbow SET® Device Options

Masimo PN	Description	UM	Product Image	Notes
9081	<p>Field Upgrader</p> <p>Used to load optional Rainbow parameters and features onto compatible monitors. Current available options include: Total Hemoglobin (SpHb), Oxygen Content (SpOC), Methemoglobin (SpMet), Carboxyhemoglobin (SpCO), Pleth Variability Index (PVI), Desat Index 3D alarm and PI Delta 3D alarm. Each optional parameter and feature sold separately.</p>	EA		<p>Maximum quantity of 25 of any individual parameter or feature per Field Upgrader.</p> <p>Includes: 1 Field Upgrader, 1 to 5 Operator's Manuals, parameter upgrade indicator sticker(s), and postage paid envelope to return Field Upgrader to Masimo.</p>
2630	<p>Total Hemoglobin (SpHb) parameter</p> <p>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field</p>	EA	No Image	
2296	<p>Carboxyhemoglobin (SpCO) parameter</p> <p>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field</p>	EA	No Image	
2297	<p>Methemoglobin (SpMet) parameter</p> <p>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field</p>	EA	No Image	
2313	<p>Pleth Variability Index (PVI)</p> <p>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field</p>	EA	No Image	
2940	<p>Acoustic Respiration (RRa)</p> <p>Can be ordered on a Rad87 (part number 9226, 9227, 9228, 9229) or on a Field Upgrader (part number 2926) to upgrade a compatible monitor in the field</p>	EA	No Image	
3518	<p>Respiration Rate from Pleth (RRp)</p> <p>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field. In countries where RRp is cleared.</p>	EA	No Image	
3519	<p>Respiration Rate from Pleth (RRp), Fixed Term License</p> <p>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field. In countries where RRp is cleared.</p>	EA	No Image	

Masimo Rainbow SET® Device Options

Masimo PN	Description	UM	Product Image	Notes
2631	Carboxyhemoglobin (SpCO), Fixed Term License <i>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field</i>	EA	No Image	
2632	Methemoglobin (SpMet), Fixed Term License <i>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field</i>	EA	No Image	
2633	Total Hemoglobin (SpHb), Fixed Term License <i>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field</i>	EA	No Image	
2636	PI Delta 3D Alarm, Fixed Term License <i>Can be ordered on a Radical-7 (part number 9086 or 9153) or on a Field Upgrader (part number 9081) to upgrade a compatible monitor in the field.</i>	EA	No Image	
2637	Desat Index 3D Alarm, Fixed Term License <i>Can be ordered on a Radical-7 (part number 9086 or 9153) or on a Field Upgrader (part number 9081) to upgrade a compatible monitor in the field.</i>	EA	No Image	
3402	Adaptive Threshold Alarm (ATA), Fixed Term License <i>Can be ordered on a Rad-87 (part numbers 9226, 9227, 9228, 9229). In countries where ATA is cleared.</i>	EA	No Image	
3325	Rainbow Acoustic Monitoring (RRa), Fixed Term License <i>Can be ordered on a Radical-7 or Rad-87 or on a Field Upgrader to upgrade a compatible monitor in the field.</i>	EA	No Image	
2762	Full Rainbow® License (SpHb, PVI, SpCO, SpMet) <i>Can be ordered on a Radical-7 (part number 9086 or 9153) or on a Field Upgrader (part number 9081) to upgrade a compatible monitor in the field.</i>	EA	No Image	
2946	Rainbow® Activation License <i>Allows Radical (20-pin) part number 9230 to be upgraded to Rainbow parameters.</i>	EA	No Image	
2341	Service Tool, MX-1	EA	No Image	Upgrades up to 100 instruments
2897	Service Tool, MX-3	EA	No Image	Upgrades up to 100 instruments
3351	Service Tool, MX-3 with ATA <i>In countries where ATA is cleared.</i>	EA	No Image	Upgrades up to 100 instruments
2340	Demonstration Tool <i>Used to exercise monitors alarm default setting and demonstrate optional parameters (SpHb, SpCO, SpMet and PVI) when optional parameters are installed on a compatible monitor.</i>	EA		

Philips® Compatible Devices

Masimo PN	Description	UM	Product Image	Notes
2267	Masimo SET® IntelliVue® Module <i>Note: Patient Cable must be ordered separately.</i>	EA		
1881	Masimo SET® VueLink® Module <i>Comes with VueLink cable.</i>	EA		
1841	VueLink® Cable <i>For use with Masimo SET VueLink Module, 3 ft.</i>	EA		




Philips®, IntelliVue® and VueLink® are registered trademarks of Philips Corporation.

SafetyNet™ System				
Masimo PN	Description	UM	Product Image	Notes
9095	SafetyNet™ System <i>Configurable for customer needs.</i>	System		
2146	SafetyNet™ Compatible Radical Docking Station, RDS-1 <i>Comes with power cord and communication capabilities.</i> <i>Note: SatShare is not supported on RadNet Compatible Radical Docking Stations.</i>	EA		
2147	SafetyNet™ Compatible Radical Docking Station, RDS-3 <i>Comes with power cord and communication capabilities.</i> <i>Note: SatShare is not supported on RadNet Compatible Radical Docking Stations.</i>	EA		
14422	3Com Workgroup Bridge <i>Wi-Fi, 802.11a/b/g Wireless Lan</i>	EA		(1) Required per Pulse Oximeter.
15007	Colubris Workgroup Bridge	EA		(1) Required per Pulse Oximeter; alternative for 3Com Workgroup Bridge.
14423	3Com Ethernet Workgroup Bridge Antenna <i>Optional accessory to 14422.</i>	EA	No Image	
15275	NET232-DTE (male) RS232 Serial to Ethernet Intelligent Cable Adapter	EA		
14530	Null Modem Cable, DB-9 Male to DB-9 Male, 1' Length	EA	No Image	

SafetyNet™ System				
Masimo PN	Description	UM	Product Image	Notes
32054	Bracket A Main (RADICAL-7) <i>Includes bracket and hardware for mounting the 3Com Radio to the Radical-7 Pulse Oximeter.</i>	EA		
32055	Bracket B Main (RAD-8) <i>Includes bracket and hardware for mounting the 3Com Radio to the Rad-8 Pulse Oximeter.</i>	EA		
14256	Cisco 1242 AP Wireless Access Point <i>Includes Cisco Access Point, power cables and basic mount. Must also include (1) 14260 Antenna for each access point.</i>	EA		
14260	Antenna, Ceiling Mount <i>2dBi DIVERSITY Wireless Antenna for 14256 Cisco Wireless Access Point (802.11 b/g).</i>	EA		(1) 14260 or 14261 Antenna required per 14256 Cisco Wireless Access Point.
14261	Antenna, Wall Mount <i>2.4GHz 6.5 dBi Wireless Antenna for 14256 (802.11 b/g).</i>	EA		May be required based on Wireless Site Survey.
14572	Cisco ANT-5135D-R Antenna <i>(802.11 a).</i>	EA		
14573	Antenna, 4.5 dBi Diversity, Cisco ANT-51 <i>Use with 14256 (802.11a Networks).</i>	EA		
14574	Antenna, 9.5 Dbi Diversity, Cisco ANT-51 <i>Use with 14256 (802.11a Networks).</i>	EA		
14858	10/100 Ethernet Switch, SMC Networks SMC6128PL2 TigerSwitch	EA		
33612	SafetyNet™ Appliance	EA		





SafetyNet™ System

Masimo PN	Description	UM	Product Image	Notes
14804	HP R1500 G2 UPS	EA		
14857	XRackPro2 Rackmount 4U Cabinet, Noise Reduction Enclosure	EA		
14868	SafetyNet™ Terminal PC	EA		
16170	SafetyNet™ View	EA		
14424	Paging Transmitter <i>Frequency Range (450 - 470 MHz).</i>	EA		(1) required per Floor in Hospital
14253	Pager 456 - 460 MHz, Synthesized, POCSAG <i>Pager tunable in the 456 - 460 MHz Range.</i>	EA		
14254	Pager 460 - 464 MHz, Synthesized, POCSAG <i>Pager tunable in the 460 - 464 MHz Range.</i>			
14255	Pager 464 - 468 MHz, Synthesized, POCSAG <i>Pager tunable in the 464 - 468 MHz Range.</i>			



SEDLine™ Brain Function Monitor				
Masimo PN	Description	UM	Product Image	Notes
9185	SEDLine Monitor	EA		
2481	SEDLine Patient Module <i>18ft. For use with SEDLine monitor.</i>	EA		
2480	SEDLine Patient Cable <i>6ft. For use with SEDLine EEG Sensors.</i>	EA		
2479	SEDLine EEG Sensor <i>Single patient use 25/box, Non-sterile This product does not contain natural rubber latex.</i>	Box		

Sensors & Cables



rainbow® Acoustic Monitoring Sensor and Cables (for use with Masimo rainbow Acoustic Monitoring Devices with RRA parameters)

Masimo PN	Description	UM	Product Image	Notes
2759	RAS-125, Acoustic Respiration Sensor Adult Adhesive Sensors (RRA) Single patient use 10/box, Non-sterile Weight > 30kg This product does not contain natural rubber latex.	Box		
2902	RAS-125c, Acoustic Respiration Cloth Sensor Adult Adhesive Sensors (RRA) Single patient use 10/box, Non-sterile Weight > 30kg This product does not contain natural rubber latex.	Box		
3475	RAS-125c, Acoustic Respiration Cloth Sensor Adult/Pediatric Adhesive Sensors (RRA) Single patient use 10/box, Non-sterile Weight > 10kg This product does not contain natural rubber latex. Available in countries where cleared.	Box		
3401	RAS-125c, Short Term Monitoring Acoustic Respiration Cloth Sensor Adult Adhesive Sensors (RRA) Single patient use 10/box, Non-sterile Weight > 30kg This product does not contain natural rubber latex. Available in the US only.	Box		
3483	RAS-125c, Short Term Monitoring Acoustic Respiration Cloth Sensor Adult/Pediatric Adhesive Sensors (RRA) Single patient use 10/box, Non-sterile Weight > 10kg This product does not contain natural rubber latex. Available in countries where cleared.	Box		





rainbow® Acoustic Monitoring Sensor and Cables (for use with Masimo rainbow Acoustic Monitoring Devices
with RRa parameters)







Masimo PN	Description	UM	Product Image	Notes
2757	Dual-Channel rainbow® Acoustic Monitoring Cable Channel 1 - Accepts any Masimo 20-pin cable or direct-connect sensor Channel 2 - Accepts the RAM Patient Cable only (PN 2758).	EA		
2758	Acoustic Respiration Patient Cable rainbow 14-pin Patient Cable, 10 ft. For use with Acoustic Respiration sensors.	EA		

rainbow® Adhesive Sensors (for use with Masimo rainbow SET Devices with SpHb, SpO₂, and SpMet parameters)





Masimo PN	Description	UM	Product Image	Notes
2414	rainbow® R1 25L Adult Adhesive Sensors (SpHb, SpO ₂ , SpMet) Single patient use 10/box, Non-sterile Weight >30 kg This product does not contain natural rubber latex. Use replacement tape 2623.	Box		Requires rainbow RC patient cable
2415	rainbow® R1 20L Infant Adhesive Sensors (SpHb, SpO ₂ , SpMet) Single patient use 10/box, Non-sterile Weight 3 kg – 30 kg This product does not contain natural rubber latex. Use replacement tape 2624.	Box		Requires rainbow RC patient cable

rainbow® Adhesive Sensors (for use with Masimo rainbow SET Devices with SpO₂, SpCO, and SpMet parameters)




Masimo PN	Description	UM	Product Image	Notes
2221	rainbow® R25 Adult Adhesive Sensors (SpO ₂ , SpCo, SpMet) Single patient use 10/box, Non-sterile Weight >30 kg This product does not contain natural rubber latex.	Box		Requires rainbow RC patient cable
2222	rainbow® R20 Pediatric Adhesive Sensors (SpO ₂ , SpCo, SpMet) Single patient use 10/box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex.	Box		Requires rainbow RC patient cable
2219	rainbow® R25-L Adult / Neonatal Adhesive Sensors (SpO ₂ , SpCo, SpMet) Single patient use 10/box, Non-sterile Weight <3 kg or >30 kg This product does not contain natural rubber latex. Use replacement tape 2623.	Box		Requires rainbow RC patient cable
2220	rainbow® R20-L Infant Adhesive Sensors (SpO ₂ , SpCo, SpMet) Single patient use 10/box, Non-sterile Weight 3 kg – 30 kg This product does not contain natural rubber latex. Use replacement tape 2624.	Box		Requires rainbow RC patient cable

rainbow® ReSposable Sensors (for use with Masimo rainbow SET Devices with SpHb, SpO ₂ , and SpMet parameters)				
Masimo PN	Description	UM	Product Image	Notes
3457	rainbow® ReSposable™ R2-25 Sensor System Adult R2-25a disposable optical sensors & R2-25r reusable optical sensors (SpO ₂ , SpHb, SpMet) Single patient use 20-R2-25a /box (DOS) 1-R2-25r /box (ROS) Non-sterile, Weight > 30 kg This product does not contain natural rubber latex.	Box		Requires rainbow RC patient cable
3458	rainbow® ReSposable™ R2-20 Sensor System Pediatric R2-20a disposable optical sensors & R2-20r reusable sensors (SpO ₂ , SpHb, SpMet) Single patient use 20-R2-20a /box (DOS) 1-R2-20r /box (ROS) Non-sterile, Weight 10 kg – 50 kg This product does not contain natural rubber latex.	Box		Requires rainbow RC patient cable
3395	rainbow® ReSposable™ R2-25r Sensors Adult R2-25r reusable optical sensors (SpO ₂ , SpHb, SpMet) Single patient use 5-R2-25r /box (ROS) Non-sterile, Weight > 30 kg This product does not contain natural rubber latex.	EA		For use with R2-25a
3411	rainbow® ReSposable™ R2-20r Sensors Pediatric/Slender Digit R2-20r reusable optical sensors (SpO ₂ , SpHb, SpMet) Single patient use 5-R2-20r /box (ROS) Non-sterile, Weight 10 – 50 kg This product does not contain natural rubber latex.	EA		For use with R2-20a
2753	rainbow® ReSposable™ R2-25a Sensors Adult R2-25a disposable optical sensors (SpO ₂ , SpHb, SpMet) Single patient use 25-R2-25a /box (DOS) Non-sterile, Weight > 30 kg This product does not contain natural rubber latex.	Box		Requires R2-25r
2755	rainbow® ReSposable™ R2-20a Sensors Pediatric/Slender Digit R2-20a disposable optical sensors (SpO ₂ , SpHb, SpMet) Single patient use 25-R2-20a /box (DOS) Non-sterile, Weight 10 – 50 kg This product does not contain natural rubber latex.	Box		Requires R2-20r



rainbow® Reusable Spot Check Sensors (for use with Masimo rainbow SET Devices with SpHb, SpO₂, SpMet parameters; SpHb parameter required)



Masimo PN	Description	UM	Product Image	Notes
3418	rainbow® DCI SC 200 Adult reusable sensor, 3 ft. (SpHb, SpO ₂ and SpMet) Weight > 30 kg This product does not contain natural rubber latex. 200 Spot Checks	EA		Requires rainbow RC Patient Cable
3419	rainbow® DCIP SC 200 Pediatric reusable sensor, 3 ft. (SpHb, SpO ₂ and SpMet) Weight 10 – 50 kg This product does not contain natural rubber latex. 200 Spot Checks	EA		Requires rainbow RC Patient Cable
3420	rainbow® DCI SC 400 Adult reusable sensor, 3 ft. (SpHb, SpO ₂ and SpMet) Weight > 30 kg This product does not contain natural rubber latex. 400 Spot Checks	EA		Requires rainbow RC Patient Cable
3421	rainbow® DCIP SC 400 Pediatric reusable sensor, 3 ft. (SpHb, SpO ₂ and SpMet) Weight 10 – 50 kg This product does not contain natural rubber latex. 400 Spot Checks	EA		Requires rainbow RC Patient Cable



Pronto-7 Reusable Spot Check Sensors

Masimo PN	Description	UM	Product Image	Notes
9400	rainbow® 4D DC Sensor, Small Reusable Sensor 1/box, Non-sterile Used with finger sizes between 8mm – 12mm diameter. This product does not contain natural rubber latex. Comes with 500 or 1000 Spot checks.	EA		Must add Spot checks separately (P/N 3312) Requires a Pronto-7.
9401	rainbow® 4D DC Sensor, Medium Reusable Sensor 1/box, Non-sterile Used with finger sizes between 12mm – 16mm diameter. This product does not contain natural rubber latex. Comes with 500 or 1000 Spot checks.	EA		Must add Spot checks separately (P/N 3312) Requires a Pronto-7.
9402	rainbow® 4D DC Sensor, Large Reusable Sensor 1/box, Non-sterile Used with finger sizes between 16mm – 20mm diameter. This product does not contain natural rubber latex. Comes with 500 or 1000 Spot checks.	EA		Must add Spot checks separately (P/N 3312) Requires a Pronto-7.








rainbow® Reusable Sensors (for use with Masimo rainbow SET Devices with SpCO, SpO₂, and SpMet parameters)

Masimo PN	Description	UM	Product Image	Notes
2696	rainbow® DCI Adult Reusable Sensor, 3 ft. (SpCO, SpMet, SpO ₂) 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		Requires rainbow RC Patient Cable
2697	rainbow® DCIP Pediatric Reusable Sensor, 3 ft. (SpCO, SpMet, SpO ₂) 1/box, Non-sterile Weight 10 kg - 50 kg This product does not contain natural rubber latex.	EA		Requires rainbow RC Patient Cable



rainbow® Direct Connect Sensors (for use with Masimo rainbow SET Devices with SpCO, SpO ₂ , and SpMet parameters)				
Masimo PN	Description	UM	Product Image	Notes
2201	rainbow® DCI-dc3 Adult Reusable Direct Connect Sensor, 3 ft. (SpO ₂ , SpCO, SpMet) 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2407	rainbow® DCI-dc8 Adult Reusable Direct Connect Sensor, 8 ft. (SpO ₂ , SpCO, SpMet) 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.			
2202	rainbow® DCI-dc12 Adult Reusable Direct Connect Sensor, 12 ft. (SpO ₂ , SpCO, SpMet) 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.			
2069	rainbow® DCIP-dc3 Pediatric Reusable Direct Connect Sensor, 3 ft. (SpO ₂ , SpCO, SpMet) 1/box, Non-sterile Weight 10 kg - 50 kg This product does not contain natural rubber latex.	EA		
2640	rainbow® DCIP-dc8 Pediatric Reusable Direct Connect Sensor, 8 ft. (SpO ₂ , SpCO, SpMet) 1/box, Non-sterile Weight 10 kg - 50 kg This product does not contain natural rubber latex.			
2070	rainbow® DCIP-dc12 Pediatric Reusable Direct Connect Sensor, 12 ft. (SpO ₂ , SpCO, SpMet) 1/box, Non-sterile Weight 10 kg - 50 kg This product does not contain natural rubber latex.			

rainbow® Patient Cables (for use with Masimo rainbow SET Devices)				
Masimo PN	Description	UM	Product Image	Notes
2405	rainbow® RC-1 <i>rainbow 20-pin Patient Cable, 1 ft</i> <i>For use with rainbow and M-LNCS sensors</i>	EA		
2406	rainbow® RC-4 <i>rainbow 20-pin Patient Cable, 4 ft</i> <i>For use with rainbow and M-LNCS sensors</i>	EA		
2404	rainbow® RC-12 <i>rainbow 20-pin Patient Cable, 12 ft</i> <i>For use with rainbow and M-LNCS sensors</i>			







Masimo SET ReSposable™ S-DOS Series

Masimo PN	Description	UM	Product Image	Notes
2434	ReSposable™ S-DOS™ 25 Adult SpO ₂ disposable optical sensors (DOS) Single patient use 20 /box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	Box		For use with S-ROS 3U or 12U series.
2435	ReSposable™ S-DOS™ 20 Pediatric SpO ₂ disposable optical sensors (DOS) Single patient use 20 /box, Non-sterile Weight 10 - 50 kg This product does not contain natural rubber latex.	Box		For use with S-ROS 3U or 12U series
2469	ReSposable™ S-DOS™ 25L Neonatal/Adult SpO ₂ disposable optical sensors (DOS) Single patient use 20 /box, Non-sterile Weight > 30 kg, < 3 kg This product does not contain natural rubber latex. Use replacement tape 3439.	Box		For use with S-ROS 3U or 12U series
2470	ReSposable™ S-DOS™ 20L Infant SpO ₂ disposable optical sensors (DOS) Single patient use 20 /box, Non-sterile Weight 3-30 kg This product does not contain natural rubber latex. Use replacement tape 3439.	Box		For use with S-ROS 3U or 12U series
2471	ReSposable™ S-DOS™ 30L Large Adult SpO ₂ disposable optical sensors (DOS) Single patient use 20 /box, Non-sterile Weight > 40 kg This product does not contain natural rubber latex. Use replacement tape 3439.	Box		For use with S-ROS 3U or 12U series
2476	ReSposable™ S-DOS™ 25Pt-L Neonatal SpO ₂ disposable optical sensors (DOS) Single patient use 20 /box, Non-sterile Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 3440.	Box		For use with S-ROS 3U or 12U series
2477	ReSposable™ S-DOS™ 20L-500 Non-adhesive Neonatal SpO ₂ disposable optical sensors (DOS) Single patient use 20 /box, Non-sterile Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 3441.	Box		For use with S-ROS 3U or 12U series




Masimo SET ReSposable™ S-ROS 3U Series







Masimo PN	Description	UM	Product Image	Notes
2473	LNCS® ReSposable™ S-ROS™ 3U (S-ROS™) <i>SpO₂ reusable optical sensor (ROS)</i> <i>Reusable</i> <i>1/box, Non-sterile</i> <i>This product does not contain natural rubber latex</i>	Box		For use with S-DOS series. Requires patient cable.
2472	M-LNCS™ ReSposable™ S-ROS™ 3U (S-ROS™) <i>SpO₂ reusable optical sensor (ROS)</i> <i>Reusable</i> <i>1/box, Non-sterile</i> <i>This product does not contain natural rubber latex</i>	Box		For use with S-DOS series. Requires patient cable.

Masimo SET ReSposable™ S-ROS 12U Series




Masimo PN	Description	UM	Product Image	Notes
2474	ReSposable™ Red S-ROS™ 12U (S-ROS™) <i>SpO₂ reusable optical sensor (ROS), 12 ft., 20 pin</i> <i>Reusable</i> <i>1/box, Non-sterile</i> <i>This product does not contain natural rubber latex</i>	Box		For use with S-DOS series.
2475	ReSposable™ S-ROS™ 12U (S-ROS™) <i>SpO₂ reusable optical sensor (ROS), 12 ft., 14 pin</i> <i>Reusable</i> <i>1/box, Non-sterile</i> <i>This product does not contain natural rubber latex</i>	Box		For use with S-DOS series.
2619	ReSposable™ SpaceLabs S-ROS™ 12U (S-ROS™) <i>SpO₂ reusable optical sensor (ROS), 12 ft.</i> <i>Reusable</i> <i>1/box, Non-sterile</i> <i>This product does not contain natural rubber latex</i>	Box		For use with S-DOS series.
2620	ReSposable™ Philips S-ROS™ 12U (S-ROS™) <i>SpO₂ reusable optical sensor (ROS), 12 ft.</i> <i>Reusable</i> <i>1/box, Non-sterile</i> <i>This product does not contain natural rubber latex</i>	Box		For use with S-DOS series.
2621	ReSposable™ DB9 S-ROS™ 12U (S-ROS™) <i>SpO₂ reusable optical sensor (ROS), 12 ft.</i> <i>Reusable</i> <i>1/box, Non-sterile</i> <i>This product does not contain natural rubber latex</i>	Box		For use with S-DOS series.
2622	ReSposable™ GE S-ROS™ 12U (S-ROS™) <i>SpO₂ reusable optical sensor (ROS), 12 ft.</i> <i>Reusable</i> <i>1/box, Non-sterile</i> <i>This product does not contain natural rubber latex</i>	Box		For use with S-DOS series.

ReSposable™ MAC S-ROS 12U Series




Masimo PN	Description	UM	Product Image	Notes
2566	ReSposable™ MAC 395 S-ROS™ 12U (S-ROS™) N-395 SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1/box, Non-sterile This product does not contain natural rubber latex	Box		For use with S-DOS series.
2567	ReSposable™ MAC 180 S-ROS™ 12U (S-ROS™) N-180 SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1/box, Non-sterile This product does not contain natural rubber latex	Box		For use with S-DOS series.
2568	ReSposable™ MAC GE S-ROS™ 12U (S-ROS™) GE SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1/box, Non-sterile This product does not contain natural rubber latex	Box		For use with S-DOS series.
2569	ReSposable™ MAC SL S-ROS™ 12U (S-ROS™) SpaceLabs SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1/box, Non-sterile This product does not contain natural rubber latex	Box		For use with S-DOS series.
2570	ReSposable™ MAC SL2 S-ROS™ 12U (S-ROS™) SpaceLabs SpO ₂ reusable optical sensor (ROS), 12 ft. Reusable 1/box, Non-sterile This product does not contain natural rubber latex	Box		For use with S-DOS series.
2571	ReSposable™ MAC CMS S-ROS™ 12U (S-ROS™) Philips SpO ₂ reusable optical sensor (ROS), 12 ft., 12 pin Reusable 1/box, Non-sterile This product does not contain natural rubber latex	Box		For use with S-DOS series.







LNCS® Adhesive Sensors				
Masimo PN	Description	UM	Product Image	Notes
1859	LNCS® Adtx Adult SpO ₂ Adhesive Sensors, 18 in. Single Patient Use 20/box, Sterilizable Weight > 30 kg This product does not contain natural rubber latex.	Box		
2317	LNCS® Adtx-3 Adult SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight > 30 kg This product does not contain natural rubber latex.			
2918	LNCS® E1™ Adult SpO ₂ disposable ear sensor, 3 ft. Single Patient Use 10/box, non-sterile, non-adhesive Weight > 30 kg This product does not contain natural rubber latex.	Box		Requires Masimo SET MS-2000 (Version 4.8 or higher) technology or Masimo rainbow SET MX technology
1860	LNCS® Pdtx Pediatric SpO ₂ Adhesive Sensors, 18 in. Single Patient Use 20/box, Sterilizable Weight 10 kg – 50 kg This product does not contain natural rubber latex.	Box		
2318	LNCS® Pdtx-3 Pediatric SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight 10 kg – 50 kg This product does not contain natural rubber latex.			
1861	LNCS® Inf-L Infant SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight 3 kg – 20 kg This product does not contain natural rubber latex.	Box		
1862	LNCS® Neo-L Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight < 3 kg or > 40 kg This product does not contain natural rubber latex.	Box		
1901	LNCS® NeoPt-L Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex.	Box		







LNCS® Adhesive Sensors with Replaceable Tapes



Masimo PN	Description	UM	Product Image	Notes
2328	LNCS® Inf Infant SpO ₂ Adhesive Sensors, 18 in. Single Patient Use 20/box, Sterilizable Weight 3 kg – 20 kg This product does not contain natural rubber latex. Use replacement tape 2307.	Box		
2319	LNCS® Inf-3 Infant SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight 3 kg – 20 kg This product does not contain natural rubber latex. Use replacement tape 2307.			
2329	LNCS® Neo Neonatal SpO ₂ Adhesive Sensors, 18 in. Single Patient Use 20/box, Sterilizable Weight < 3 kg or > 40 kg This product does not contain natural rubber latex. Use replacement tape 2308.	Box		
2320	LNCS® Neo-3 Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight < 3 kg or > 40 kg This product does not contain natural rubber latex. Use replacement tape 2308.			
2330	LNCS® NeoPt Neonatal SpO ₂ Adhesive Sensors, 18 in. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2321	LNCS® NeoPt-3 Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2309.			
2331	LNCS® NeoPt-500 Neonatal SpO ₂ non-adhesive sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2322.	Box		




LNCS® Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat setting





Masimo PN	Description	UM	Product Image	Notes
2411	LNCS® Trauma Adult SpO ₂ Adhesive Sensors Single Patient Use 20/box, Sterilizable Weight > 30 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2412	LNCS® Newborn Neonatal Neonatal SpO ₂ Adhesive Sensors Single Patient Use 20/box, Sterilizable Weight < 3 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2413	LNCS® Newborn Infant/Pediatrics Infant / Pediatric SpO ₂ Adhesive Sensor Single Patient Use 20/box, Sterilizable Weight 3 kg – 30 kg This product does not contain natural rubber latex. Use replacement wrap 2322.	Box		

LNCS® Sterile Sensors				
Masimo PN	Description	UM	Product Image	Notes
2601	LNCS® Aadx-3 Adult SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterile Weight > 30 kg This product does not contain natural rubber latex.	Box		Note: The sensors are <u>not</u> intended to be used for evaluations or carried as sales samples. If there is a need for samples of sterile sensors, contact your manager for approval and include the e-mail with FRM-1386.
2602	LNCS® Pdtx-3 Pediatric SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex.	Box		
2603	LNCS® Inf-3 Infant SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterile Weight 3 kg – 20 kg This product does not contain natural rubber latex. Use replacement tape 2307.	Box		
2604	LNCS® Neo-3 Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterile Weight < 3 kg or > 40 kg This product does not contain natural rubber latex. Use replacement tape 2308.	Box		
2605	LNCS® NeoPt-3 Neonatal SpO ₂ Adhesive Sensors, 3 ft. Single Patient Use 20/box, Sterile Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2606	LNCS® NeoPt-500 Neonatal SpO ₂ non-adhesive sensor, 18 in. Single Patient Use 20/box, Sterile Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2322.	Box		

LNCS® Reusable Sensors				
Masimo PN	Description	UM	Product Image	Notes
1863	LNCS® DCI Adult Reusable Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
1864	LNCS® DCIP Pediatric Reusable Sensor, 3 ft. 1/box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex.	EA		
2653	LNCS® DBI Adult Soft Reusable Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
1895	LNCS® TC-I Reusable Tip-Clip EAR Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2258	LNCS® YI Multisite Reusable Sensor, 3 ft. 1/box, Non-sterile Multiple Foam and Adhesive Wraps Weight > 1 kg This product does not contain natural rubber latex.	EA		
1896	LNCS® TF-I Adult Transflectance Forehead, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		

LNCS® Red Patient Cables (for use with Masimo Rainbow SET Devices)				
Masimo PN	Description	UM	Product Image	Notes
2365	Red LNC-01 <i>LNCS 20-pin SpO₂ Spot Check Patient Cable.</i>	EA		
2055	Red LNC-04 <i>LNCS 20-pin SpO₂ Patient Cable, 4 ft.</i>	EA		
2056	Red LNC-10 <i>LNCS 20-pin SpO₂ Patient Cable, 10 ft.</i>			
2057	Red LNC-14 <i>LNCS 20-pin SpO₂ Patient Cable, 14 ft.</i>			

LNC™ Patient Cables (for use with Masimo SET Devices)				
Masimo PN	Description	UM	Product Image	Notes
2364	LNC-1 <i>LNCS Series Patient Cable.</i>	EA		
2017	LNC-4 <i>LNCS Series Patient Cable, 4 ft.</i>	EA		
1814	LNC-10 <i>LNCS Series Patient Cable, 10 ft.</i>			
2013	LNC-14 <i>LNCS Series Patient Cable, 14 ft.</i>			
2021	LNC-4-Ext <i>LNC Extension Cable, 4 ft.</i>	EA		
2281	LNC MP10 <i>LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO2 LNCS Dual Key Patient Cable, 10 ft.</i>	EA		

LNC™ Patient Cables (for use with Masimo SET Devices)				
Masimo PN	Description	UM	Product Image	Notes
1902	LNC-10 P8 <i>LNCS to IntelliVue Masimo SET Patient Cable, 10 ft.</i>	EA		
2430	LNC-4-GE <i>LNCS to GE Patient Cable, 4ft.</i>	EA		
2016	LNC-10-GE <i>LNCS to GE Patient Cable, 10ft.</i>			
2432	LNC-SL-10 <i>LNCS to Spacelabs Patient Cable, 10 ft.</i>	EA		
2285	LNC-04-DB9 <i>LNCS to DB9 Connector Patient Cable, 4 ft.</i>	EA		
2286	LNC-10-DB9 <i>LNCS to DB9 Connector Patient Cable, 10 ft.</i>			
2287	LNC-14-DB9 <i>LNCS to DB9 Connector Patient Cable, 14 ft.</i>			




LNCS® Adapter Cables

Masimo PN	Description	UM	Product Image	Notes
1816	LNCS® to PC Adapter Cable <i>LNCS Series to LNOP PC Adapter Cable, 18 in.</i>	EA		
2263	LNC™ CMS <i>LNCS Series to Philips 12-pin Instrument Adapter Cable, 3 m.</i>	EA		
2264	LNC™ MAC-GE <i>LNCS Series to GE Adapter Cable, 3 m.</i>	EA		
2266	LNC™ MAC-SL <i>LNCS Series to Spacelabs Adapter Cable, 3 m.</i>	EA		
2336	LNC™ MAC-SL2 <i>LNCS Series to Spacelabs Adapter Cable, 3 m.</i>	EA		

LNCS® Adapter Cables

Masimo PN	Description	UM	Product Image	Notes
2268	LNC™ MAC-180 <i>LNCS Series to N-180 Adapter Cable, 10 ft.</i>	EA		
2269	LNC™ MAC-395 <i>LNCS Series to N-395 Adapter Cable, 10 ft.</i>	EA		
2921	LNC-MAC-4 P8 <i>LNCS Series to Philips FAST-SpO2 Adapter Cable, 4 ft.</i>	EA		
2922	LNC-MAC-10 P8 <i>LNCS Series to Philips FAST-SpO2 Adapter Cable, 10 ft.</i>			


M-LNCS™ Adhesive Sensors







Masimo PN	Description	UM	Product Image	Notes
2508	M-LNCS™ Adtx Adult SpO ₂ Adhesive Sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight > 30 kg This product does not contain natural rubber latex.	Box		
2509	M-LNCS™ Adtx-3 Adult SpO ₂ Adhesive Sensor, 3 ft. Single Patient Use 20/box, Sterilizable Weight > 30 kg This product does not contain natural rubber latex.			
2919	M-LNCS™ E1™ Adult SpO ₂ disposable ear sensor, 3 ft. Single Patient Use 10/box, non-sterile, non-adhesive Weight > 30 kg This product does not contain natural rubber latex.	Box		Requires Masimo SET MS-2000 (Version 4.8 or higher) technology or Masimo rainbow SET MX technology
2510	M-LNCS™ Pdtx Pediatric SpO ₂ Adhesive Sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight 10 kg – 50 kg This product does not contain natural rubber latex.	Box		
2511	M-LNCS™ Pdtx-3 Pediatric SpO ₂ Adhesive Sensor, 3 ft. Single Patient Use 20/box, Sterilizable Weight 10 kg – 50 kg This product does not contain natural rubber latex.			

M-LNCS™ Adhesive Sensors with Replaceable Tapes




Masimo PN	Description	UM	Product Image	Notes
2512	M-LNCS™ Inf Infant SpO ₂ Adhesive Sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight 3 kg – 20 kg This product does not contain natural rubber latex. Use replacement tape 2307.	Box		
2513	M-LNCS™ Inf-3 Infant SpO ₂ Adhesive Sensor, 3 ft. Single Patient Use 20/box, Sterilizable Weight 3 kg – 20 kg This product does not contain natural rubber latex. Use replacement tape 2307.			
2514	M-LNCS™ Neo Neonatal/Adult SpO ₂ Adhesive Sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight < 3 kg or > 40 kg This product does not contain natural rubber latex. Use replacement tape 2308.	Box		
2515	M-LNCS™ Neo-3 Neonatal/Adult SpO ₂ Adhesive Sensor, 3 ft. Single Patient Use 20/box, Sterilizable Weight < 3 kg or > 40 kg This product does not contain natural rubber latex. Use replacement tape 2308.			
2516	M-LNCS™ NeoPt Neonatal SpO ₂ Adhesive Sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2517	M-LNCS™ NeoPt-3 Neonatal SpO ₂ Adhesive Sensor, 3 ft. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2309.			
2518	M-LNCS™ NeoPt-500 Neonatal SpO ₂ Non-Adhesive Sensor, 18 in. Single Patient Use 20/box, Sterilizable Weight < 1 kg This product does not contain natural rubber latex. Use replacement wrap 2322.	Box		

M-LNCS™ Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat setting.



Masimo PN	Description	UM	Product Image	Notes
2521	M-LNCS™ Trauma Adult SpO ₂ Sensor Single Patient Use 20/box, Sterilizable Weight > 30 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2519	M-LNCS™ Newborn Neonatal Neonatal SpO ₂ Sensor Single Patient Use 20/box, Sterilizable Weight < 3 kg This product does not contain natural rubber latex. Use replacement wrap 2309.	Box		
2520	M-LNCS™ Newborn Infant/Pediatric Infant/Pediatric SpO ₂ Sensor Single Patient Use 20/box, Sterilizable Weight 3 kg – 30 kg This product does not contain natural rubber latex. Use replacement wrap 2322.	Box		
2522	M-LNCS™ Blue Neonatal/Infant/Pediatric SpO ₂ Adhesive Sensor Single Patient Use 10/box, Non-Sterile Weight 2.5 – 30kg This product does not contain natural rubber latex. Use replacement tape 2822.	Box		

M-LNCS™ Reusable Sensors				
Masimo PN	Description	UM	Product Image	Notes
2501	M-LNCS™ DCI Adult Reusable Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2502	M-LNCS™ DCIP Pediatric/Slender Digit Reusable Sensor, 3 ft. 1/box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex.	EA		
2507	M-LNCS® DBI Adult Soft Reusable Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2503	M-LNCS™ TC-I Reusable Tip-Clip EAR Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2505	M-LNCS™ YI Multisite Reusable Sensor, 3 ft. 1/box, Non-sterile Multiple Foam and Adhesive Wraps Weight > 1 kg This product does not contain natural rubber latex.	EA		
2504	M-LNCS™ TF-I Adult Transflectance Forehead, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		






M-LNC™ Patient Cables

Masimo PN	Description	UM	Product Image	Notes
2523	M-LNC-1 <i>M-LNCS Series Patient Cable, 1 ft.</i>	Box		
2524	M-LNC-4 <i>M-LNCS Series Patient Cable, 4 ft.</i>	Box		
2525	M-LNC-10 <i>M-LNCS Series Patient Cable, 10 ft.</i>			
2526	M-LNC-14 <i>M-LNCS Series Patient Cable, 14 ft.</i>			
2831	M-LNC-MP10 <i>M-LNCS to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO2 M-LNC Dual Key Patient Cable, 10 ft.</i>	Box		
2532	M-LNC-10 P8 <i>M-LNCS to IntelliVue Masimo SET M-LNC Patient Cable, 10 ft.</i>	Box		
2655	M-LNC-4-GE <i>M-LNCS to GE Patient Cable, 4 ft.</i>	Box		
2533	M-LNC-10-GE <i>M-LNCS to GE Patient Cable, 10 ft.</i>			

M-LNC™ Patient Cables

Masimo PN	Description	UM	Product Image	Notes
2714	M-LNC-SL-10 <i>M-LNCS to Spacelabs Patient Cable, 10 ft.</i>	Box		
2534	M-LNC-04-DB9 <i>M-LNCS Series to DB9 Connector Patient Cable, 4 ft.</i>	Box		
2535	M-LNC-10-DB9 <i>M-LNCS Series to DB9 Connector Patient Cable, 10 ft.</i>			
2536	M-LNC-14-DB9 <i>M-LNCS Series to DB9 Connector Patient Cable, 14 ft.</i>			






M-LNCS™ Adapter Cables

Masimo PN	Description	UM	Product Image	Notes
2542	M-LNC™ MAC-SL2 <i>M-LNC Series to Spacelabs Adapter Cable, 2 ft.</i>	EA		
2541	M-LNC™ MAC-SL <i>M-LNC Series to Spacelabs Adapter Cable, 3 m.</i>	EA		
2537	M-LNC™ MAC-180 <i>M-LNCS Series to N-180 Adapter Cable, 10 ft.</i>	EA		
2538	M-LNC™ MAC-395 <i>M-LNCS Series to N-395 Adapter Cable, 10 ft.</i>	EA		
2748	M-LNC™ MAC P8 <i>M-LNCS to Philips FAST-SpO2 Adapter Cable, 10 ft.</i>	EA		





M-LNCS™ Adapter Cables

Masimo PN	Description	UM	Product Image	Notes
2540	M-LNC™ MAC-GE <i>M-LNCS Series to GE Adapter Cable, 3 m.</i>	EA		
2539	M-LNC™ MAC CMS <i>M-LNCS Series to Phillips 12-pin Instruments Adapter Cable, 3 m.</i>	EA		
2654	M-LNCS™ to LNC Adapter Cable <i>M-LNCS Series to LNC Patient Cable, 1.5 ft.</i>	EA		
2750	M-LNCS™ to LNC-MAC Adapter Cable <i>M-LNCS Series to LNC-MAC Adapter Cable, 1.5 ft.</i>	EA		

LNOP® Adhesive Sensors

Masimo PN	Description	UM	Product Image	Notes
1001	LNOP® Adt Adult SpO ₂ Adhesive Sensors Single Patient Use 20/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	Box		
1025	LNOP® Pdt Pediatric SpO ₂ Adhesive Sensors Single Patient Use 20/box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex.	Box		
1800	LNOP® Inf-L Infant SpO ₂ Adhesive Sensors Single Patient Use, 20/box, Non-sterile Weight 3 kg – 20 kg This product does not contain natural rubber latex. Use 1801 replacement tape.	Box		
1798	LNOP® Neo-L Neonatal SpO ₂ Adhesive Sensors Single Patient Use, 20/box, Non-sterile Weight < 3 kg or > 40 kg This product does not contain natural rubber latex. Use 1799 replacement tape.	Box		
1651	LNOP® NeoPt-L Neonatal SpO ₂ Adhesive Sensors Single Patient Use 20/box, Non-sterile Weight < 1 kg This product does not contain natural rubber latex. Use 1926 replacement tape.	Box		

LNOP® Specialty Sensors enables the algorithms to give new data higher priority than old data regardless of the averaging, sensitivity, or FastSat setting.

Masimo PN	Description	UM	Product Image	Notes
2358	LNOP® Trauma Adult SpO ₂ Adhesive Sensors Single Patient Use 20/box, Non-sterile Weight >30 kg This product does not contain natural rubber latex. Use 1926 replacement tape	Box		
2203	LNOP® Newborn Neonatal Neonatal SpO ₂ Adhesive Sensors Single Patient Use 20/box, Non-sterile Weight < 3 kg This product does not contain natural rubber latex. Use 1926 replacement tape	Box		
2204	LNOP® Newborn Infant/ Pediatric Infant/ Pediatric SpO ₂ Adhesive Sensors Single Patient Use, 20/box, Non-sterile Weight 3 kg – 30 kg This product does not contain natural rubber latex. Use 2391 replacement tape	Box		
1970	LNOP® Blue Neonatal /Infant/ Pediatric SpO ₂ Adhesive Sensors Single Patient Use 10/box, Non-sterile Weight 2.5 kg – 30 kg This product does not contain natural rubber latex. Use 1801 replacement tape.	Box		

LNOP® Reusable Sensors







Masimo PN	Description	UM	Product Image	Notes
1269	LNOP® DCI Adult SpO ₂ Reusable Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
1276	LNOP® DCIP Pediatric SpO ₂ Reusable Sensor, 3 ft. 1/box, Non-sterile Weight 10 kg – 50 kg This product does not contain natural rubber latex.	EA		
1794	LNOP® TC-I SpO ₂ Reusable Tip-Clip EAR Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
1793	LNOP® TF-I Adult SpO ₂ Transflectance Forehead Reusable Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
1560	LNOP® DC-195 Adult SpO ₂ Reusable Sensor 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
1544	LNOP® YI Multisite SpO ₂ Reusable Sensor, 3 ft. 1/box, Non-sterile Multiple Foam and Adhesive Wraps Weight > 1 kg This product does not contain natural rubber latex.	EA		

LNOP® Red Patient Cables (for use with Masimo Rainbow SET Devices)



Masimo PN	Description	UM	Product Image	Notes
2058	Red PC-04 LNOP Series 20-pin SpO ₂ Patient Cable, 4 ft.	EA		
2059	Red PC-08 LNOP Series 20-pin SpO ₂ Patient Cable, 8 ft.			
2060	Red PC-12 LNOP Series 20-pin SpO ₂ Patient Cable, 12 ft.			

LNOP® Patient Cables (for use with Masimo SET Devices)				
Masimo PN	Description	UM	Product Image	Notes
1173	PC04 LNOP Series Patient Cable, 4 ft.	EA		
1005	PC08 LNOP Series Patient Cable, 8 ft.			
1006	PC12 LNOP Series Patient Cable, 12 ft.			
1619	PC04-Ext LNOP Extension Cable, 4 ft.	EA		
2282	LNOP® MP12 LNOP to IntelliVue Masimo SET or IntelliVue Philips FAST-SpO2 LNOP Dual Key Patient Cable, 12 ft.	EA		
1847	LNOP® PC-12 P8 LNOP to Masimo SET Patient Cable, 12 ft.	EA		
1890	PC-12-GE LNOP to GE Patient Cable, 12ft.	EA		
2288	PC-04-DB9 LNOP Series to DB9 Connector Patient Cable, 4ft.	EA		
2289	PC-08-DB9 LNOP Series to DB9 Connector Patient Cable, 10ft.			
2290	PC-12-DB9 LNOP Series to DB9 Connector Patient Cable, 12 ft.			
2433	LNOP® PC-SL-12 LNOP Spacelabs Patient Cable, 12ft.	EA		




LNOP® Adapter Cables

Masimo PN	Description	UM	Product Image	Notes
1848	LNOP® MAC P8 <i>LNOP to Philips FAST-SpO2 Adapter Cable, 10 ft.</i>	EA		
1849	LNOP® CMS <i>LNOP Series to Philips 12-pin Instruments Adapter Cable, 12 ft.</i>	EA		
1645	MAC-1 <i>LNOP Series to LNC Adapter Cable, 18 in.</i>	EA		
2858	M-LNCS™ to PC Adapter Cable <i>M-LNCS Series to LNOP PC Patient Cables, 18 in.</i>	EA		
2014	Drager Patient Cable <i>For use with Masimo SET LNOP sensors, 2m</i>	EA		
1927	MAC-1 <i>LNOP SpO2 Patient Cable, 12 ft.</i>	EA		





Direct Connect Sensors (for use with Masimo SET Devices)

Masimo PN	Description	UM	Product Image	Notes
1396	DCSC Adult Reusable Direct Connect SpO ₂ Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2387	DC-8 Adult Reusable Direct Connect SpO ₂ Sensor, 8 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.			
1969	DC-12 Adult Reusable Direct Connect SpO ₂ Sensor, 12 ft. 1/box, Non-sterile. Weight > 30 kg This product does not contain natural rubber latex.			
2651	DBI-dc3 Adult Soft Reusable Direct Connect SpO ₂ Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2652	DBI-dc8 Adult Soft Reusable Direct Connect SpO ₂ Sensor, 8 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.			








Red™ Direct Connect Sensors (for use with Masimo Rainbow SET Devices)

Masimo PN	Description	UM	Product Image	Notes
2053	Red DCI-dc3 Adult Reusable Direct Connect SpO ₂ Sensor, 3 ft. 1/box, Non-sterile Weight > 30kg This product does not contain natural rubber latex.	EA		
2054	Red DCI-dc12 Adult Reusable Direct Connect SpO ₂ Sensor, 12 ft. 1/box, Non-sterile Weight > 30kg This product does not contain natural rubber latex.			
2256	Red DCIP-dc3 Pediatric Reusable Direct Connect SpO ₂ Sensor, 3 ft. 1/box, Non-sterile Weight 10kg – 50kg This product does not contain natural rubber latex.	EA		
2257	Red DCIP-dc12 Pediatric Reusable Direct Connect SpO ₂ Sensor, 12 ft. 1/box, Non-sterile Weight 10kg – 50kg This product does not contain natural rubber latex.			
2643	Red DBI-dc3 Adult Soft Reusable Direct Connect Sensor, 3 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.	EA		
2644	Red DBI-dc8 Adult Soft Reusable Direct Connect Sensor, 8 ft. 1/box, Non-sterile Weight > 30 kg This product does not contain natural rubber latex.			








Mountable Patient Cable Extensions and Adapter

Masimo PN	Description	UM	Product Image	Notes
1520	PXC16 Mountable Patient Cable Extension, 16 ft. For use with Masimo SET LNC, M-LNC or PC Series Patient Cables.	EA		Discontinued. Replacement p/n's: 2371 and 2373
1555	PXC30 Mountable Patient Cable Extension, 30 ft. For use with Masimo SET LNC, M-LNC or PC Series Patient Cables.			
1557	PXC50 Mountable Patient Cable Extension, 50 ft. For use with Masimo SET LNC, M-LNC or PC Series Patient Cables.			
2478	14 pin to 20 pin Adapter Cable Connect PXC2 to Masimo SET® 20 pin devices, SpO2 only.	EA		
2371	PXC2-16 Mountable Extension Cable, 16-ft. For use with Masimo SET® PC, LNC™, M-LNC™ Patient Cables and DCSC, DC-8, DC-12, DBI-dc3, DBI-dc8 Direct Connect sensors or Masimo 14 pin to 20 pin Adapter Cable.	EA		
2373	PXC2-30 Mountable Extension Cable, 30-ft. For use with Masimo SET® PC, LNC™, M-LNC™ Patient Cables and DCSC, DC-8, DC-12, DBI-dc3, DBI-dc8 Direct Connect sensors or Masimo 14 pin to 20 pin Adapter Cable.	EA		









Accessories

Sensor Accessories				
Masimo PN	Description	UM	Product Image	Notes
2624	Replacement Tapes Rainbow® R1 20L and R20-L For use with Rainbow R1 20L and R20-L 50/box, Non-sterile This product does not contain natural rubber latex.	EA		
2623	Replacement Tapes For Rainbow® R1 25L and R25-L For use with Rainbow R1 25L and R25-L 50/box, Non-sterile This product does not contain natural rubber latex.	EA		
1799	Replacement Tapes For LNOP® Neo-L Series Sensors For use with LNOP Neo-L 100/pack, Non-sterile This product does not contain natural rubber latex.	EA		
1801	Replacement Tapes For LNOP® Inf-L and LNOP Blue Sensors For use with LNOP Inf-L and LNOP Blue 100/pack, Non-sterile This product does not contain natural rubber latex.	EA		
1926	Replacement Wraps for LNOP® NeoPt-L, LNOP® Trauma and Newborn Neonatal Sensors For use with LNOP NeoPt-L, Trauma and Newborn Neonatal 10/pack, Non-sterile This product does not contain natural rubber latex.	EA		
2391	Replacement Wraps for LNOP® Newborn Infant / Pediatric Sensors For use with LNOP Newborn Infant / Pediatric 10/pack, Non-sterile This product does not contain natural rubber latex.	EA		
2307	Replacement Tapes for LNCS®/M-LNCS™ Inf Series Sensors For use with LNCS/M-LNCS Inf/Inf-3 102/box, Non-sterile This product does not contain natural rubber latex.	EA		

Sensor Accessories




Masimo PN	Description	UM	Product Image	Notes
2308	Replacement Tapes For LNCS®/M-LNCS™ Neo Series Sensors For use with LNCS/M-LNCS Neo/Neo-3 102/box, Non-sterile This product does not contain natural rubber latex.	EA		
2309	Replacement Wraps for LNCS®/M-LNCS™ NeoPt, NeoPt-3, Trauma and Newborn Neonatal Sensors For use with LNCS/M-LNCS NeoPt/ NeoPt-3/ Trauma and Newborn Neonatal 10/pack, Non-sterile This product does not contain natural rubber latex.	EA		
2322	Replacement SofTouch Wraps for LNCS®/M-LNCS™ NeoPt-500, Newborn Infant / Pediatric For use with LNCS/M-LNCS NeoPt-500/Newborn Infant/Pediatric 10/pack, Non-sterile This product does not contain natural rubber latex.	EA		
2822	Replacement Tape for M-LNCS™ Blue For use with M-LNCS Blue 50/pack, Non-sterile This product does not contain natural rubber latex.	EA		
3439	Replacement Tape for ReSposable S-DOS 30L, S-DOS 25L and S-DOS 20L For use with ReSposable S-DOS 30L/S-DOS 25L/S-DOS 20L 100/box, Non-sterile This product does not contain natural rubber latex.	EA		
3440	Replacement Tape for ReSposable S-DOS 25Pt-L For use with ReSposable S-DOS 25Pt-L 10/pack, Non-sterile This product does not contain natural rubber latex.	EA		
3441	Replacement Tape for ReSposable S-DOS 20L-500 For use with ReSposable S-DOS 20L-500 10/pack, Non-sterile This product does not contain natural rubber latex.	EA		

Sensor Accessories

Masimo PN	Description	UM	Product Image	Notes
1596	CleanShield® Multisite Wrap For use with YI Sensors Adult/Pediatric/Neonatal 100/box, Non-sterile This product does not contain natural rubber latex.	EA		
1597	Standard Multisite Wrap For use with YI Sensors Adult/Pediatric/Neonatal 100/box, non-sterile This product does not contain natural rubber latex.	EA		
1598	Standard Petite Wrap For use with YI Sensors Adult slender digit/Ped/Neonatal adhesive 100/box, non-sterile This product does not contain natural rubber latex.	EA		
1600	Accessory Adhesive Squares For use with YI Sensors 144/pack, non-sterile	EA		
1903	Transflectance Sensor Adhesive Pad For use with TF-I Sensors 20/pack, non-sterile	EA		
16179	Tele Plast Philips SpO2 Connector Plug. Compatible with Philips MMS module.	EA		
1602	Foam Wraps for LNOP®/LNCS®/M-LNCS™ YI Sensors For use with YI Sensors 12/pack, non-sterile Washable This product does not contain natural rubber latex.	EA		
1601	Accessory Clothing Clips For securing reusable sensor cables to patient clothing 5/pack, non-sterile	EA		

Sensor Accessories










Masimo PN	Description	UM	Product Image	Notes
3444	Masimo® Soft Wrap Armband Used to secure the sensor cable along the patients arm. 20/pack, non-sterile This product does not contain natural rubber latex.	EA		
1608	Masimo Headband For use with TF-I Sensors 5/box, non-sterile This product does not contain natural rubber latex.	EA		
2215	Masimo Headband For use with TF-I Sensors 12/box, non-sterile This product does not contain natural rubber latex.			
21514	LNC™ Replacement Cover For use with all LNC patient cables 25 covers and 50 pins.	EA		
3443	Masimo® Cable Tether Secures reusable sensor to patient cable. 20/pack, non-sterile This product does not contain natural rubber latex.	EA		




Training Kits				
Masimo PN	Description	UM	Product Image	Notes
3454	SpO₂ ReSpposable, Adult/Pediatric Sensor Training Kit <i>Comes with 1 S-DOS 25, 1 S-DOS 20, 1 S-DOS 25L, 1 S-DOS 20L, 1 S-DOS 30L, 1 soft wrap armband, 1 cable tether, and 1 sensor application poster, direction for use.</i>	EA		
3455	SpO₂ ReSpposable, L Series Sensor Training Kit <i>Comes with 1 S-DOS 30L, 1 S-DOS 25L, 1 S-DOS 20L, 1 S-DOS 25Pt-L, 1 soft wrap armband, 1 cable tether, 1 sensor application poster and 1 neonatal foot with application card, direction for use.</i>	EA		
2659	M-LNCS® Sensor Training Kit <i>Includes 1 LNCS Amtx, 1 LNCS Pdtx, 1 LNCS Inf, 1 LNCS Neo, 1 LNCS NeoPt sensor and 1 neonatal foot with application card, direction for use.</i>	EA		

Training Kits






Masimo PN	Description	UM	Product Image	Notes
2345	LNCS® Sensor Training Kit <i>Includes 1 LNCS Adtx, 1 LNCS Pdtx, 1 LNCS Inf, 1 LNCS Neo, 1 LNCS NeoPt sensor and 1 neonatal foot with application card, direction for use.</i>	EA		
1282	LNOP® Sensors Application Training Packet <i>Comes with 1 LNOP Adt, 1 LNOP Pdt, 1 LNOP Inf-L, 1 LNOP Neo-L and 1 LNOP NeoPt-L, 2 cotton swabs, 2 alcohol swabs and 1 neonatal foot with application card, direction for use.</i>	EA		
3456	Acoustic Respiration Sample Pack <i>1-RAS-125/1-RAS-125C Sensor</i>	EA		

Bedside Device Accessories				
Masimo PN	Description	UM	Product Image	Notes
1315	Radical® Series Replacement Battery <i>For use with Radical and Radical-7Blue Screen Handhelds.</i>	EA		
2386	Radical-7™ Color Screen Replacement Battery <i>For use with Radical-7 Color Screen only.</i>	EA		
1317	Radical® Series Pole Clamp <i>For use with Radical and Radical-7 Handhelds.</i>	EA		
1395	Radical® Series Handheld Lock <i>To secure Radical and Radical-7 Blue Screen Handhelds in Radical Docking Stations.</i>	EA		
30534	Radical® Series Handheld Lock Key <i>To unlock Radical and Radical-7 Blue Screen handhelds from Docking Stations.</i>	EA		
1584	Radical® Series Power Cord Lock <i>To secure a power cord to Radical Docking Stations 5/pack.</i>	EA		
2228	Rad-8® Power Cord Lock	EA		
30475	Radical® Series Handheld Case <i>(Shown with Blue Screen Handheld Radical. Blue Radical not included).</i>	EA		








Bedside Device Accessories				
Masimo PN	Description	UM	Product Image	Notes
2367	9-pin Masimo SET® Tester <i>Tester for devices containing Masimo SET and connects to LNCS patient cables (Tests only SpO2 and BPM).</i>	EA		
2368	Masimo Rainbow SET® Tester <i>Tester for devices containing Masimo SET or Masimo Rainbow SET, 20 pin instrument connector. (Tests SpHb, SpCO, SpMet, PVI, PI, SpO2 and BPM).</i>	EA		Can be issued with all 20 pin devices
1593	14 Pin Masimo SET® Tester <i>For Oximeter SpO2 and PR Testing Tester for devices containing Masimo SET, 14 pin instrument connector. (Tests only SpO2 and BPM).</i>	EA		
1795	Masimo SET® Tester <i>For Oximeter & Patient Cable SpO2 and PR Testing. Tester for devices containing Masimo SET and connects to LNOP PC Patient Cables. (Tests only SpO2 and BPM).</i>	EA		
2227	Rad-8® Mounting Bracket	EA		
2350	Radical-7 Color Handheld Lock <i>Locks Radical-7 Color handheld to the docking station (used with PN 2351).</i>	EA		
2351	Radical-7 Color Handheld Lock Key <i>Used for installation of Radical-7 Color Handheld Lock and unlocking handheld from docking station (Used with PN 2350)</i>	EA		
2389	Mini-Demo Connector - 30 Min <i>Used to demo product (best used to display trending information).</i>	EA		
2390	Mini-Demo Connector - 2 Min <i>Used to demo product (best used for short demonstrations).</i>	EA		

Handheld Device Accessories				
Masimo PN	Description	UM	Product Image	Notes
2207	Water Resistant Protective Carrying Case for Battery operated portable handhelds, Black <i>(This can be used with Rad-5/5v).</i>	EA	 <p><i>Shown with Quick Reference Card, Rainbow® Sensor and Rad-57. Not included.</i></p>	
2208	Water Resistant Protective Carrying Case for Battery operated portable handhelds, Red <i>(This can be used with Rad-5/5v).</i>	EA	 <p><i>Shown with Quick Reference Card, Rainbow® Sensor and Rad-57. Not included.</i></p>	
13158	Handheld Carrying Case with Shoulder Strap <i>(Shown with Rad-57. Device not included).</i>	EA		Note: Compatible with all Masimo Handheld Devices (Rad-5, Rad-5v, Rad-57, Pronto).

Handheld Device Accessories

Masimo PN	Description	UM	Product Image	Notes
10368	Radical-7 Ambulation Pouch <i>Pack of 10</i>	EA	No Image	
1909	Data Transfer Download Cable <i>Used for connecting Rad-5 (with 14-pin connector) to Personal Computer to enable download of trending data.</i>	EA		
2063	Data Transfer Download Cable, Red – 20 Pin <i>Used for connecting Rad-57's and 20-pin Rad-5's to Personal Computer to enable download of trending data.</i>	EA		
1908	TrendCom™, Trend Download Software <i>Used with appropriate Data Transfer Download cable to initiate download and de-compression of Rad-5 and Rad-57 trending data to a Personal Computer. Used with serial cable to initiate download and decompression of Radical-7, Rad-87 and Rad-8 trending data.</i>	EA		
1842	Handheld Boot, Gray <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		
1980	Handheld Boot, Yellow <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		

Handheld Device Accessories

Masimo PN	Description	UM	Product Image	Notes
1981	Handheld Boot, Red <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		
1982	Handheld Boot, Orange <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		
2097	Handheld Boot, Royal Blue <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		
2098	Handheld Boot, Light Blue <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		
2099	Handheld Boot, Pink <i>For use with Rad-5, Rad-5v and all Rad-57 handhelds.</i>	EA		
2357	Light Shield <i>For use in strobing environment and high sunlight intensity areas. 5/Pack</i>	EA		
2815	Masimo Ambient Shield Accessory <i>For use with Rainbow Sensors. 10/pack</i>	EA		











Pronto-7™ Spot Check Pulse CO-Oximeter Accessories

Masimo PN	Description	UM	Product Image	Notes
3382	Boot, Pronto-7 Red	EA		
3383	Boot, Pronto-7 Yellow	EA		
3384	Boot, Pronto-7 Green	EA		
3385	Boot, Pronto-7 Light Blue	EA		
3386	Boot, Pronto-7 Pink	EA		








Pronto-7™ Spot Check Pulse CO-Oximeter Accessories

Masimo PN	Description	UM	Product Image	Notes
3387	Boot, Pronto-7 Dark Blue	EA		
3388	Boot, Pronto-7 Black	EA		
3389	Boot, Pronto-7 Orange	EA		
34983	Finger Gauge, Pronto-7	EA		
16782	Travel Bag, Pronto-7	EA		
16868	Bracket, Pronto-7	EA		








SatShare® Cables

Masimo PN	Description	Validated Monitor Interfaces [Plug-In Parameter Modules]*	UM	Host Connector	Notes
1595	AT01	Airshields Athena Monitor with Nellcor SpO ₂	EA		
1324	DO01	Datex-Ohmeda AS/3 Patient Monitor [M-NSAT] System 5 Compact Monitor [M-NSAT]	EA		
1326	DO02	Datex-Ohmeda AS/3 Patient Monitor [M-NESTPR-00-00 / M-NESTR-00-00] CardioCap/5 Monitor Light Monitor System 5 Compact Monitor [M-NESTPR-00-00 / M-NESTR-00-00]	EA		
1528	DO03	Datex-Ohmeda AS/3 Patient Monitor [M-EST]	EA		
1533	DO04	Datex-Ohmeda AS/3 Patient Monitor [M-ESTP] System 5 Compact Monitor [M-EST / M-ESTP]	EA		
2023	DO05	Datex-Ohmeda System 5 Compact Monitor [M-NESTPR-00-01]	EA		
2223	DO06	Datex-Ohmeda System 5 Compact Monitor [M-PRESTN..01]	EA		
2969	DO07	Datex-Ohmeda System 5 Anesthesia Monitor [E-PRESTN]	EA		
1325	DS01	Corometrics Maternal / Fetal Monitor Model 118 Maternal / Fetal Monitor Model 120 Critikon Dinamap MPS Select Monitor Datascope Expert Monitor with Nellcor SpO ₂ [HS-500] Accutor Monitor with Nellcor SpO ₂ Passport Monitor with Nellcor SpO ₂ North American Draeger Evita 2 Dura, Evita 4 Anesthesia Machine Julian Anesthesia Machine Narcomed 4 Anesthesia Machine Vitalert 1000, 3100 Monitor Edentec / Nellcor Puritan Bennett Edentrace II Plus (Model 3711) Medical Data Electronics Escort Prism Monitor [20001] Escort / Escort II Monitor Welch Allyn Protocol Propaq CS Monitor Propaq Encore Monitor Propaq Series 100 Monitor Witt Biomedical Frontend Patient Monitor	EA		
1539	DS02	Datascope Passport Monitor with Datascope SpO ₂ Accutor Monitor with Datascope SpO ₂	EA		

SatShare® Cables

Masimo PN	Description	Validated Monitor Interfaces [Plug-In Parameter Modules]*	UM	Host Connector	Notes
1365	DS03	Datascope Expert Fukuda Denshi Dynascope DS-5300W [HS-500]	EA		
2979	DS04	Datascope Fukuda Denshi Dynascope DS-5400 [HS-500]	EA		
1357	HP03	Philips / Agilent / HP Viridia CMS Model [M1020A] M1165, M1166, M1167, M1175, M1176, M1177 Monitors [M1020A] M1092, M1094, M1095, M1097 Display Modules, using the M1046 Computer Module [M1020A] Omnicare / Viridia Model 24/26 Series Monitor System M1204 Model 78352 Monitor [M1020A] Series 50XM Fetal Monitor	EA		
1321	MQ01	GE / Marquette Solar 7000 / 8000 Series Monitor, Tramscope 12 / 12C Monitor, Series 7000 / 7000RA / 7005 / 7010 / 7010RA / 7200 Monitor, Smart- Pac Transport Display Monitor [Tram 250SL / Tram 450SL / Tram 650SL / Tram 850A&SL / SpO2 Solar Parameter Module] Eagle 3000 Series Monitor Eagle 4000 Series Monitor Dash Series Monitor North American Draeger Narkomed 6100 Anesthesia Systems with Integrated Patient Monitoring	EA		
1789	MS01	Datex-Ohmeda S5 [E-Masimo-00] Draeger Infinity Delta [SmartPOD with Masimo SET] GE/Marquette Solar 8000 [Tram 451M / Tram 851M / Masimo SET single parameter module] Philips IntelliVue Masimo SET MMS SpO2 module	EA		
1367	NK01	Nihon Kohden Model BMS-8301A/MV-831RA and Model BMS-8502A [AL-800PA]	EA		
1368	NK02	Nihon Kohden BSM-2301A and BSM-5105K	EA		

SatShare® Cables

Masimo PN	Description	Validated Monitor Interfaces [Plug-In Parameter Modules]*	UM	Host Connector	Notes
1366	OH01	Airshields Athena Monitor with Ohmeda SpO ₂ Ohmeda RGM Monitor (Respiratory Gas Machine) GE / Marquette Solar 7000 / 8000 Series Monitor, Tramscope 12 / 12C Monitor, Series 7000 / 7000RA / 7005 / 7010 / 7010RA / 7200 Monitor, Smart-Pac Transport Display Monitor [Tram 100A / Tram 200A&SL / Tram 400A / Tram 500A / Tram 600A&SL / Tram 800A&SL]	EA		
1646	OH02	Hellige SMU 611 Monitor <i>(Not available for sale in the United States)</i>	EA		
1840	PM01	Philips Medical M3/M4 [MMS Module] Intellivue Monitor [MMS Module]	EA		
1323	SL01	SpaceLabs Medical PC Express Monitor (90308), PC Scout Monitor (90309), PC Bedside / Central Monitor (90303, 90311, 90312), PC2 Bedside / Central Monitor (90305, 90315, 90316) [90465 / 90466 / 90467 / 90489] 600 Series Monitor	EA		
1362	SL02	SpaceLabs Medical Ultraview 1030 / 1050 Portable Bedside Monitor (90367, 90369) [90496] Ultraview 1500 / 1600 / 1700 Bedside Monitor (90363, 90364, 90387) [90496] Ultraview UCW Bedside Monitor (90385) [90496]	EA		
1322	SM01	Siemens Medical Systems Sirecust 1281 Monitor SC 6002XL, SC7000, SC9000XL	EA		
1605	SatShare Extension Cable		EA		

* Module compatibility information is subject to change without notice. For more current information, please contact Masimo Corporation.